

# EXO GAG

## SERUM

Next generation of Exosome isolation  
for biomarker discovery in liquid biopsy

### EXOSOME PURIFICATION KIT FOR SERUM

#### USER GUIDE

#### STORAGE

All components can be stored at room temperature.

#### PRODUCT COMPONENT

EXOGAG - Serum™ isolation kit 50ml (50-25 samples). 1 x User guide.  
Not supplied: collection tubes.

**EXOGAG - SERUM is a specific, quick and inexpensive method to optimize the process of exosome isolation.**

#### PRODUCT INFORMATION

The EXOGAG - Serum precipitation reaction is based on the interaction between the precipitation solution and glycosaminoglycans located on exosome surface. After a simple incubation, the exosomes can be isolated by a short centrifugation.


#### WHAT IS EXOGAG - SERUM

EXOGAG - Serum technology allows the isolation of exosomes from serum based on the affinity of exosomes to the precipitation reagent, which allows their isolation from a complex sample. EXOGAG - Serum is a patented exosome purification method that allows the isolation of exosomes from a serum sample with a minimal amount of co-precipitated material, such as protein or genetic material (DNA, RNA or microRNA) which makes it an ideal product for biomarkers research and their transfer to the clinic.



High sensibility and specificity	Fast	Inexpensive	Small sample needed	No specific equipment needed
----------------------------------	------	-------------	---------------------	------------------------------

SAMPLE	SAMPLE, vol.	EXOGEN - SERUM, vol.
Serum	500µl	1ml
Serum	1ml	2ml



## PROTOCOL

### EXOGEN - Serum Exosomes Isolation Protocol.

- 1. Collect serum sample.** Samples can be frozen until the moment they are used; if the samples have been frozen, thaw and temper them before processing.
- 2. Centrifuge the sample** at 2000 x g for 5' to remove cells and cell debris.
- 3. Transfer the supernatant** to a new tube and discard the pellet of possible cell debris.
- 4. Add the volume of sample** to isolate exosomes to a new tube and add twice the volume of EXOGEN - Serum precipitation reagent, as shown in the table.
- 5. Mix the sample** and EXOGEN - Serum precipitation reagent by inverting the tube or vortexing to homogenize the final solution (the solution will have a characteristic blue colour).
- 6. Incubate the sample** for 5' at 4°C.
- 7. Centrifuge the sample** at 3000 x g; 15' at 4°C.
- 8. Remove the supernatant** being careful not to remove the pellet containing the exosomes (this pellet will be dark blue).
- 9. Resuspend the exosomes** in the appropriate buffer (pipetting repeatedly up and down), depending on the technique.